

# Optical Microscopes

# Lexicon


❖ **Micro:** very small

microgram, micrometer, microsecond, micron, microdot, microbiology, microorganism, microsurgery, microeconomics, microelectronics, microwave, microchip, microcomputer, microprocessor, microphone, microfilm, .....

❖ **Microscope:** an instrument for magnifying very small objects


❖ **Microscopic:** so small as to be visible only with a microscope


❖ **Microscopy:** the use of a microscope



# Microscope

One or more lens that makes an enlarged image of an object.





Oldest published image known to have been made with a microscope: Bees by **Francesco Stelluti**, 1630.

- ❖ Optical microscope, often referred to as “**light microscope**”, is a type of microscope which **uses visible light** & a system of lenses to magnify images of small samples.
- ❖ Optical microscopes are the oldest design of microscope & were possibly designed in their present compound form in 17<sup>th</sup> century.
- ❖ Aim: to improve **RESOLUTION** & **CONTRAST**.
- ❖ Microscopes which do not use visible light are:
  - Scanning Electron Microscope (SEM)**
  - Transmission Electron Microscope (TEM)**
  - Atomic Force Microscope (AFM)**

- ❖ Inventor of microscope: **Galileo Galilei**

Galileo developed a compound microscope with a convex & concave lens in 1609.

- ❖ **Giovanni Faber** coined the name “microscope”.

Greek words:

**micron** meaning “small”



**skopein** meaning “to look at”

- ❖ Optical & electron microscopy involve **diffraction**, **reflection**, or **refraction** of electromagnetic radiation/electron beams interacting with specimen, & subsequent collection of this scattered radiation or another signal in order to create an image.
- ❖ This process may be carried out by wide-field irradiation of sample (e.g., standard **light microscopy** & **transmission electron microscopy**) or by scanning of a fine beam over sample (e.g., **confocal laser scanning microscopy** & **scanning electron microscopy**).
- ❖ **Scanning probe microscopy** involves interaction of a scanning probe with surface of object of interest.



# Lenses & Bending of Light

- ❖ Lenses focus light rays at a specific place, called focal point.
- ❖ Strength of lens is related to focal length.  
Short focal length → more magnification
- ❖ Light is refracted (bent) when passing from one medium to another.
- ❖ **Refractive index**: a measure of how greatly a substance slows velocity of light.
- ❖ Direction & magnitude of bending is determined by refractive indices between the two media forming the interface.



# Eyepiece Lens

Usually has a power of 10 X.



**Eyepiece Lens × Objective Lens = Total Magnification**

**Objective Lens:**    Low power = 4x  
                              Medium power = 10x  
                              High power = 40x

# Microscope Resolution

- ❖ Ability of a lens to separate or distinguish small objects that are close together.
- ❖ Wavelength of light used is major factor in resolution

**Shorter wavelength → Greater resolution**



## Properties of Microscope Objectives

Property	Objective			
	Scanning	Low Power	High Power	Oil Immersion
Magnification	4×	10×	40–45×	90–100×
Numerical aperture	0.10	0.25	0.55–0.65	1.25–1.4
Approximate focal length ( <i>f</i> )	40 mm	16 mm	4 mm	1.8–2.0 mm
Working distance	17–20 mm	4–8 mm	0.5–0.7 mm	0.1 mm
Approximate resolving power with light of 450 nm (blue light)	2.3 $\mu\text{m}$	0.9 $\mu\text{m}$	0.35 $\mu\text{m}$	0.18 $\mu\text{m}$


Objective

Working distance

Slide with  
specimen



## Objective Lens



# Types of Microscope

- ❖ Simple microscope
- ❖ Compound microscope
- ❖ Stereoscopic microscope
- ❖ Electron microscope
- ❖ Phase-Contrast microscope
- ❖ Digital holographic microscope

# Simple Microscope

Similar to magnifying glass & has only one lens.






# Compound Microscope

Lets light pass through an object & then through two or more lenses.



# Stereoscopic Microscope

**Gives a three-dimensional view of an object.**  
(Ex. Insects & leaves).



# Electron Microscope


- ❖ Uses a magnetic field to bend beams of electrons; instead of using lenses to bend beams of light.
- ❖ Wavelength of electron beam is much shorter than light, resulting in much higher resolution.



Range of light  
microscope

Range of  
electron microscope


Scanning  
tunneling  
microscope



# Scanning Electron Microscope

- ❖ SEM uses electrons reflected from surface of a specimen to create image.
- ❖ Sample is scanned with a beam of electrons in a raster scan pattern.
- ❖ Electrons interact with atoms that make up the sample producing signals that contain information about sample's surface topography, composition, & other properties such as electrical conductivity.


# SEM



# Transmission Electron Microscope

- ❖ Electrons scatter when they pass through thin sections of a specimen.
- ❖ Transmitted electrons (those that do not scatter) are used to produce image.
- ❖ Denser regions in specimen, scatter more electrons & appear darker.

## Transmission electron microscope







# Phase-Contrast Microscope

- ❖ Enhances contrast between intracellular structures having slight differences in refractive index.
- ❖ Excellent way to observe living cells.

## Phase-contrast microscope



## Phase-contrast microscope



# Bright-Field Microscope

- ❑ Produces a dark image against a brighter background.
- ❑ It uses several objective lenses – parfocal microscopes remain in focus when objectives are changed.
- ❑ Total magnification:  
product of magnifications of ocular lens  
& objective lens

# Dark-Field Microscope

- ❑ Produces a bright image of object against a dark background.
- ❑ It is used to observe living, unstained preparations.

# Digital Holography

# Digital Holographic Microscope

- ❑ Holography was invented by Dennis Gabor to improve electron microscope.
- ❑ Basic concept of DHM is to magnify hologram image by adopting an optical lens system so that microscope fringes can be resolved.
- ❑ DHM, unlike other microscopy, doesn't record projected image of object, rather light wavefront information originating from object is digitally recorded as a hologram.
- ❑ Imaging lens in traditional microscopy is replaced by a computer algorithm.

# Applications of DHM

**DHM has capability of non-invasively visualizing & quantifying biological tissues.**

## **Biomedical applications of DHM:**

- ❑ To perform cell counting & to measure cell viability directly in cell culture chamber.
- ❑ To study apoptotic process (programmed cell death) in different cell types. Refractive index changes taking place during apoptotic process are easily measured with DHM.
- ❑ **Cell cycle analysis:** Phase shift induced by cells has been shown to be correlated to cell dry mass, which can be combined with other parameters obtainable by DH, such as, cell volume & refractive index, to provide a better understanding of cell cycle.




- ❑ **Morphology analysis of cells:** to study cell morphology using neither staining nor labeling.
- ❑ DHM is used for automated **plant stem cell** monitoring.
- ❑ To study undisturbed processes in **nerve cells** as no labeling is required. Swelling & shape changing of nerve cells caused by cellular imbalance is easily studied.
- ❑ To measure **3-D motion of human red blood cells** moving in a microtube flow. Phase shift images are used to study **red blood cell** dynamics.
- ❑ Red blood cell volume & hemoglobin concentration are measured by combining information from absorption & phase shift images to facilitate **complete blood cell count**.
- ❑ By combining several images calculated from same hologram, but at different focal planes, an **increased depth of field** is obtained.

# Advantages

- ❑ **Simplicity of microscope:** It requires a laser, a pinhole, & a CCD camera, but no lenses at all (no aberration correction required).
- ❑ **Simplicity of sample preparation in biology:** no sectioning or staining are required, so that living cells can be viewed.
- ❑ **Maximum information:** a single hologram contains all information about 3-D structure of object.
- ❑ **Speed:** changes in specimen can ultimately be followed at capture video rate of CCD chip.
- ❑ **Maximum resolution** of order of  $\lambda$ h of laser can easily be obtained, & can be further improved by at least a factor of two or three with setup of immersion holography.
- ❑ Compared to OCT, DHM requires only a pair of particle hologram images to get complete 3D flow information.

# Digital Holography



BE: beam expander, BSs: beam splitters, SLM: spatial light modulator, RPM: random phase mask, CCD: charge coupled device, L: lens



**Phase Holographic Imaging's The Holomonitor™ M3 (Sweden)**  
[www.phiab.se](http://www.phiab.se)



**Resolutions Optics's Desktop System (Canada)**

[www.resolutionoptics.com](http://www.resolutionoptics.com)

**Submersible system** is a product with all functionality of 3D imaging technology **encased in a waterproof housing**. It allows quickly & easily observation of micro-organisms & particles up to a **depth of 5 kilometers**.




**Resolutions Optics's Submersible System (Canada)**

[www.resolutionoptics.com](http://www.resolutionoptics.com)




**Digital Holographic Microscope DHMT1000 [Lyncee tec, Switzerland]**  
[www.lynceetec.com](http://www.lynceetec.com)




**Phase-shifting digital holography.** BE: beam expander, BS: beam splitter, RP: retardation plate, M: mirror, CCD: charge-coupled device.




## Digital Hologram of Onion Peel (10X)




# Intensity of Numerical Reconstruction with DH Onion Peel (10X)



## Phase of Numerical Reconstruction with DH Onion (10X)




# 3-D presentation of Numerical Reconstruction's Phase with DH of Onion (10X)







## Digital Hologram of Onion Peel (20X)



# Intensity of Numerical Reconstruction with DH of Onion Peel (20X)




## Phase of Numerical Reconstruction with DH of Onion (20X)






## Digital Hologram of *E.coli* (20X)







# Intensity of Numerical Reconstruction with DH of *E. coli* (20X)




# Phase of Numerical Reconstruction with DH of *E. coli* (20X)



## Digital Hologram of *E.coli* (40X)



# Intensity of Numerical Reconstruction with DH of *E. coli* (40X)





## Phase of Numerical Reconstruction with DH of *E. coli* (40X)

