


**PH 301**


**ENGINEERING OPTICS**

**Lecture\_Interferometers\_16**

# What happens as light moves through scope?




# What happens as light moves through the scope?




**Microscope light is white light,  
i.e. it's made up of lots of different wavelengths;  
Each wavelength of light corresponds to a different color**


**Can prove this with a prism,  
which separates white light into its  
constituent wavelengths/colors**



# What happens as light moves through scope?




**Polarized Light**



light vibrates in all planes that contain the light ray (i.e., all planes perpendicular to the propagation direction)


**Unpolarized Light**

# 1) Light passes through the **lower polarizer**



Only the component of light vibrating in E-W direction can pass through lower polarizer – light intensity decreases


## 2) Insert the upper polarizer



Now what happens?  
What reaches your eye?

Why would anyone design a microscope that prevents light from reaching your eye?


### 3) Now insert a thin section of a rock




**How does this work?**

# Interference & Polarization


In phase




linearly polarized



Phase lag




circularly polarized






# Polarized Light Microscope

Imaging a normal sample



# Polarized Light Microscope



## Imaging a birefringent sample



- Birefringent sample splits light into e- & o- rays, which see different refractive indices
- Phase retardation of one ray with respect to other gives rise to elliptically polarized light, which is transmitted by the polarizer

# Polarized Light Microscope



Imaging a birefringent sample




Birefringent sample is bright on dark background

# Polarized Light Microscope

Add a compensator (wave plate) for better contrast




# Polarization Interferometer



# Polarizing Microscope

**Orthoscope:** An instrument used to examine the eye that eliminates corneal refraction by means of a layer of water.



**Conoscope:** A polarizing microscope for giving interference fringes & for determining principal axis of a crystal.

# Dual-Polarizing Interferometer

Dual-polarization interferometry (DPI) probes molecular layers adsorbed to the surface of a waveguide using evanescent wave of a laser beam.

It is used to measure conformational change in proteins, or other biomolecules, as they function.

DPI focuses laser light into two waveguides. One of these functions as "sensing" waveguide having an exposed surface while second one functions to maintain a reference beam.

DPI technique rotates polarization of laser, to alternately excite two polarization modes of waveguides. Measurement of interferogram for both polarizations allows both r.i. & thickness of adsorbed layer to be calculated.

Polarization can be switched rapidly, allowing real-time measurements of chemical reactions taking place on a chip surface in a flow-through system.


It is quantitative & real-time (10 Hz) with a dimensional resolution of 0.01 nm.

# Polarized Light Microscopy


- **Good for**
  - **Seeing ordered structures in the cell:**
  - **Spindles**
  - **Other cytoskeletal structures**
  - **Membranes**
  - **Collagen**
- **No staining required!**



# PBS-TGI




# Encryption using FRT & digital holography




BE: beam expander, BSs: beam splitters, SLM: spatial light modulator, RPM: random phase mask, CCD: charge-coupled device, L: lens

# Fourier domain encoding




Original image



Encrypted hologram



Encrypted image




Key hologram




Retrieval with correct keys


# Fractional Fourier domain encoding




Original image



Encrypted hologram




Encrypted image



Key hologram



Retrieval with correct keys




Retrieval with wrong fractional order



# Fully phase encryption

## Fourier domain encoding




Input image




Decrypted image




Phase-encoded image



Encrypted hologram





Key hologram



Encrypted image

# Fresnel domain encoding



Encrypted hologram



Key hologram




Encrypted image






Decrypted image


# Flexible optical encryption with multiple users & multiple security levels




# Contd...



Multiplexed encrypted image   Decrypted image - User 1   Decrypted image - User 2




Decrypted image - User n




Decryption with one of the  
wrong keys – Any User




# Experimental Set-up




# Snapshots




# Optical Recording & Reconstruction




Dice



Coin




Pin-cover




Bolt


# Optical Recording & Reconstruction




Reference Hologram



Target Hologram




Reconstructed Reference




Reconstructed Target

# Matching Holograms



**Auto-correlation**



**Cross-correlation**

**Normalized CC obtained when reference hologram is matched with target hologram using conventional JFRTC**